



User's Manual

Table of Contents

Table of Contents	1
1. Home Screen	4
1.1. Work Orders Tab	4
1.2. Sample Testing Tab	5
2. Analysis Management	6
2.2. Manage Schedules Tab	8
2.3. Analysis Scheduling Tab	9
2.4. Severity Classes Tab	11
3. Configuration Screen	13
3.1. Configuration Directory Tab	13
3.2. Connectivity Tab	14
3.3. Plants Tab	14
3.4. Groups Tab	15
3.5. Users Tab	15
3.5.1. Basic Sub-Tab	16
3.5.2. Permissions Sub-Tab	17
3.5.3. Groups Sub-Tab	19
3.5.4. User Plants Sub-Tab	19
3.6. Printers Tab	20
3.7. Settings Tab	20
3.8. Locations Tab	21
3.9. Key Shortcuts Tab	23
4. Data Management Screen	23
5. Investigation Screen	24
5.1. Configure Investigations	26
5.1.1. Constraints Sub-Tab	26
5.1.2. Alerts and Notifications Sub-Tab	26
6. Locations Screen	27
7. Product Management Screen	29
8. Report Viewer	29
8.1. Report Selector	30
8.2. Crystal Reports Viewer	30
9. Sampling Screen	31

9.1. Visualization	33
9.2. Graphing	36
9.3. Export Data	38
10. Work Orders	38
10.1. Edit Work Order Tab	38
10.2. Sample Schedule	40
11. How to insert the iPod in the Linea Pro scanner.	44
11.1. Install the ITrack Mobile App on the iPod.....	44
11.2. Perform Sample.....	47
11.3. Incubate Sample.....	48
11.4. Record Analysis.....	48
12. How to:	49
How to create a work order.....	49
How to hide unnecessary columns from the work order screen.	50
How to open a closed work order.	50
How to edit a closed sample.....	50
How to invalidate a sample.....	50
How to send an analysis for investigation.	50
How to set up a sample.....	52
An analysis disappeared from Analyses screen, and I did not delete it.	53
How to set up a threshold.	53
How to make thresholds plant-specific for global analyses.	54
How to set up a schedule.....	55
How to set up a severity/proximity/zone matrix.	56
How to add a plant.	58
How to delete a plant.....	58
How to add a user.	58
How to delete a user.....	59
How to change a user password.	59
How to create groups.	60
How to set permissions to groups.	61
How to assign users to groups.....	61
How to set permission for each user.....	62
User's user permission is different from his/her group permission.....	64

How to assign users to plants.	64
How to configure printer settings.	65
How to print labels.....	65
How to print a report.	66
How to narrow down history search using option filters.	66
How to add/edit/delete locations.	67
How to create a plant map.	69
Why aren't my locations listed in the Location drop-down list in the Work Order screen?.....	70
Why aren't all my testable locations not showing in the Location drop-down list in the Work Order screen?.....	71
How to add a product.	71
How to delete a product.	72
How to inactivate a product.	72
How to assign a product to a plant.	72
How to attach an image to a product.	72
How to apply updates.	73
My update finished with errors. What do I do?	75
How to send samples from the Sample Schedule to the work order.	76
If I add a sample to a work order manually, will it affect the status of the same sample that's in the schedule?	76
How to view/graph historical data.	76
Why isn't an analysis option showing up on my Sampling Graph?	78
How to view sample acceptability status over the plant map.	78
Part of data not being sent to plant map.....	80
How to pull up a sample on the scanner without scanning a barcode.....	80
How to create a work order on the scanner.	81
What do the result status colors mean on the scanner?	82
How to change the date range(reduce clutter) in the work order queue.	82
Recently created work order is not showing in the <i>Home Screen</i> queue.....	83
Part of data not exported to the .csv file.	84

1. Home Screen

The *Main Menu* is the command center. All of the different screens for **Presage** are available in the *Main Menu*, located on the upper left side of the screen.

Below the *Main Menu* is *Open Windows*. Every time a new screen is opened from the *Main Menu*, it will show up under the *Open Windows*. With this list, it is easy to toggle back and forth between opened screens without losing changes.

Users can easily switch from one open window to another by clicking on it or by using the **Ctrl** and **Tab** keys on the keyboard.

The top box where the Presage logo is a separate browser. There is an option to set a page with often-viewed graphs for easy monitoring. There is also an option to set any other browser page in here.

The bottom box is the work order queue. The panel on the left is a way to filter work orders, both open and closed.



1.1. Work Orders Tab

The *Work Order* tab on the *Home Screen* lists all pending work orders. Work orders are scheduled to a specific plant, date, and user. This list also specifies the date the work order was created and the unique work order number. There is also a place for notes and instructions next to each work order.



There is a menu to the left of the display panel in the *Work Orders* tab. *User* is the login user.

Date Range drop-down list shows time periods during which work orders were assigned. For example, selecting “Last Month” will list all work orders created in the previous month. There is an option to set the date range to *Custom*. For this option, users can set a date range by specifying the start and end date. To the right of the date box is a calendar icon for easier date selection.

Select the plant from the drop-down to narrow down the work order list to one plant.

By checking the *Show Closed* box, users can view all completed and pending work orders within the date range set.

Once a work order is created, open the *Home Screen* and click *Refresh* to see the new work order.

On top panel in the *Home Screen*, click *File* then *Logout* once you are finished using the program. User is automatically logged out once the program is closed or if *Exit* is selected from the *File* drop-down.

1.2. Sample Testing Tab

Sample Testing tab is next to the *Work Orders* tab in the *Home Screen*. Its purpose and function is very similar to that of the *Work Order* tab. This screen lists the testing samples which have to be pulled out of the testing mode. For example, it lists when samples have to be pulled out of the incubator.

No actions have to be taken in this screen. Once the samples are scanned to indicate when they have been pulled out, these samples are taken off of this list.

Items in this list are color-coded:

Dark grey if its been tested, Black if testing has not begun yet, Red if it should have been finished already but isn't, Blue if it should be finished 'soon' (currently in the next 24 hours), otherwise Green.


2. Analysis Management

2.1. Analyses Tab

To create a new analysis, click on the box *<Enter Analysis Name>*. Enter the analysis name on the box provided. In the next column, type in the category name. This column is used to group analyses.

Testing Time is for lab use. Type in the number hours the sample must rest in the incubator.

Keep the analysis name selected before going to the next section, *Analysis options/questions*.



Under *Option*, enter what is being tested for; for example: plate count, sample area, detect-ability, etc.

When appropriate, enter the unit in the next column; for example: cfu, ppm, etc.


The **Informational** box is for options which don't need to be graphed. Check if appropriate.

Then choose the *Value Type* from the drop-down list provided. For example, plate count is measured in number, and qualitative analyses are measured in boolean: positive or negative. To create a drop-down list of various options to select from, choose *Choice*.

Keep analysis name and analysis option selected to go on to the next section: Option choices/boundaries. Set thresholds for each option separately.

First, select the severity from the drop-down list. To set up severity classes, refer to *Analysis Management: Severity Classes*.

For *Threshold*, select an option from the drop-down list. Marginal values are those which are acceptable, but are nearing the unacceptable boundary. Unacceptable values represent failure. Invalid values are those which have no value or validity.



Constraint Type sets the minimum and maximum boundary value. Choose the appropriate *Choice/Boundary Value*, or fill it in manually. For example, to set 10 cfu as a maximum acceptable limit for an analysis, choose *Unacceptable when* under *Thresholds*; *Above* under *Constraint Type*; and type in *10* under *Choice/Boundary*.


Click **Apply** to save changes.

2.2. Manage Schedules Tab

This is where schedule intervals are defined. Choose a name for the schedule that will best identify the interval (i.e. Mondays, Weekly, By Week's End, etc.). Make sure to give each schedule a unique name that best identifies it. Click on *<Enter Schedule Name>* for the cursor to appear; type in the new schedule title.

Enter a start date, usually today's date. End Date is optional. If the end date is not set, the schedule will continue indefinitely.

Select the schedule type from the drop-down list under Type. For custom schedule, specify the number of days within one period interval in the next column. In the final column, specify the time of day (optional) the test is scheduled for.




To delete a schedule, click on the corresponding red box with white X button on the far right of each line. Click **Apply** to save changes.

2.3. Analysis Scheduling Tab

Use this screen to build an analysis schedule for each plant. First select the plant from the drop-down list.

Check the box next to *Show expired schedules* to see all expired, in addition to current, schedules.



Under the column title *Analysis*, select the analysis type from the drop-down list. This list includes all analyses created under *Analysis Management: Analyses*.

Select the schedule period interval from the drop-down list under *Schedule*. This list is created under *Analysis Management: Manage Schedules*.

The *Location Count* is automatically filled out when the locations are selected in the next step.

Analysis/Schedule:

Analysis	Schedule	Loc Count	
TPC	Semi-Monthly Pre...	103	
Enterobacteriaceae	Semi-Monthly Pre...	1	
Sanitation 1	Sanitation Pre-Op	1	
		0	



The right side of the screen is for choosing locations to perform the analyses on. Select one analysis under *Analysis/Schedule*:. Then check the box next to the locations which need to be sampled.

For a quick look-up, type in the first letter or number of the display code for the location. The corresponding locations will show up in the display below. The list can also be filtered by severity class. Choose the severity class from the drop-down; the list in the display will narrow down to locations with the severity class specified.

Add/Remove locations quick lookup:

Severity Class Filter

All



Use the **Check...** button to select groups of locations. Use the **Uncheck...** button to deselect groups of locations.

Click **Add** (on the middle of the screen, not the one on the lower portion of the screen) to add the selected locations to the list of locations to be sampled. Click **Remove** to take out locations from the list.

Code	Location	Description
<input type="checkbox"/> ST1	ST1	Stage 1
<input checked="" type="checkbox"/> 1	ST1-1	Stage 1-1
<input checked="" type="checkbox"/> 2	ST1-2	Stage 1-2
<input checked="" type="checkbox"/> 3	ST1-3	Stage 1-3
<input checked="" type="checkbox"/> 4	ST1-4	Stage 1-4
<input type="checkbox"/> ST2	ST2	Stage 2
<input checked="" type="checkbox"/> 5	ST2-5	Stage 2-5
<input checked="" type="checkbox"/> 6	ST2-6	Stage 2-6

Locations to be sampled:


Location	Plant
ST1-2	LNK - LINCOLN PLANT
ST1-3	LNK - LINCOLN PLANT
ST1-4	LNK - LINCOLN PLANT
ST8-94	LNK - LINCOLN PLANT
ST8-95	LNK - LINCOLN PLANT
ST8-96	LNK - LINCOLN PLANT
ST8-97	LNK - LINCOLN PLANT
ST8-98	LNK - LINCOLN PLANT
ST8-99	LNK - LINCOLN PLANT

To add locations one at a time, use the **Add** button on the lower portion of the screen under *Locations to be sampled*:. To remove a location, select a location from the list and click **Remove**. To remove multiple locations, hold down the *Ctrl* button on keyboard while selecting multiple locations, then click **Remove**.

Click **Apply** to save any changes and use the **Refresh/Reset** button to remove any unsaved changes.

2.4. Severity Classes Tab

This screen is for assigning severity classes to zone/proximity combinations.



First give a unique names to all severity classes. Severity classes are based on limits which are set based on product proximities and process zones. Examples are High, Medium, and Low; or 10, 100, 500. Use terminology or numbers which best fit your system.

Click on the box under *Default* (may have to expand the column to see the full column title), to set a default. Default severity class is usually the one used the most often.

Product proximity refers to contact or non-contact surfaces. To add new product proximities, click on <New Proximity>. Type in the new proximity title. To remove product proximity, click on the red box with a white X to the right of each line. To rename, select the appropriate row and type in the new proximity name.

Process zones are different operational zones such as Chilling or Packaging. To add new zones, click on <New Zone> (It may be necessary to scroll down to see this row.) To delete a zone, click on the red box with white X to the right of each line. To rename, select the appropriate row and type in the new zone name.

To assign severity classes to zone/proximity combinations, click on one severity class under *Severity Classes* on the left display to highlight it, then check the boxes where the selected severity class is appropriate on the right display.

Process Zone	Contact	Non-Contact	Environmental
Zone 1	<input checked="" type="checkbox"/> High	<input type="checkbox"/> Default	<input type="checkbox"/> Low
Zone 2	<input type="checkbox"/> Medium	<input type="checkbox"/> Default	<input type="checkbox"/> Default
Zone 3	<input checked="" type="checkbox"/> High	<input type="checkbox"/> Default	<input type="checkbox"/> Low
Zone 4	<input type="checkbox"/> Default	<input type="checkbox"/> Default	<input type="checkbox"/> Default
Zone 5	<input checked="" type="checkbox"/> High	<input type="checkbox"/> Medium	<input type="checkbox"/> Default
Zone 6	<input type="checkbox"/> Default	<input type="checkbox"/> Default	<input type="checkbox"/> Low
Zone 7	<input checked="" type="checkbox"/> High	<input type="checkbox"/> Medium	<input type="checkbox"/> Default
Zone 8	<input type="checkbox"/> Medium	<input type="checkbox"/> Default	<input type="checkbox"/> Low
Zone 9	<input checked="" type="checkbox"/> High	<input type="checkbox"/> Medium	<input type="checkbox"/> Default



Now, click the next severity class to highlight it, then check the boxes where this severity class is appropriate on the right display. Not all boxes have to be checked.


Click **Apply** to save the changes.

3. Configuration Screen

The Configuration Screen is used to configure a variety of different functions within Presage.

3.1. Configuration Directory Tab

This screen is used as a reference. The buttons here are shortcuts to the rest of the tabs in the *Configuration Screen*.





3.2. Connectivity Tab

This tab will usually be filled out by Presage Support. When appropriate, fill out the information in the space provided and click **Apply** to save the changes.

3.3. Plants Tab

Use *Plants* tab to add and delete plants.

Presage - Data Analysis System

File Edit Actions Configure Help

Main Menu

- Analysis Management
- Configuration
- Data Management
- Graphing
- Home
- Investigation
- ITrack Online
- Locations
- Report Viewer
- Sampling
- Work Orders

Open Windows

- Configuration

Configuration Directory Connectivity Graphing Configuration **Plants** Groups Users Printers Settings Locations Key Shortcuts

Plant # **New** New Plant Delete Plant

☐ Private plant (users have to log into this plant to access it)

Display Code: LNK Company Name: LINCOLN PLANT

Address: 2124 Y ST.

City: LINCOLN State: NE Zip Code: 68503 Country: USA

Phone #: Fax #

Apply Reset Changes

Click on a screen to open

Updater: (9:38 AM) software up to date | Current User: Support @ 1 NUM



To add a plant, click **New Plant**. For *Plant #*, select *New* from the drop-down list. If the new plant is private, check the box provided. With this option, users must log in to have access to the new plant.


The *Display Code* will appear in the *Plant #* drop-down list after the new plant has been created. Choose a unique code which will best identify the plant.

Company Name refers to the business name a particular plant is associated with. Fill in all information in the spaces provided.

To delete a plant, choose the plant from the *Plant #* drop-down list then click **Delete Plant**.

3.4. Groups Tab

Use Groups tab to add, delete, and set permissions to groups.



To add a group, click **Add Group**. Type in group name in the pop-up box, then click **OK**. The new group will show up under *Groups list*.


To delete a group, select the group name, then click **Delete Group**.

To set permissions to a group, select the group name, then select the permission setting from the drop-down list under *Grant* next to each category. If one permission level is appropriate for multiple categories, then select those categories while holding down the *Ctrl* key and select permission setting from the drop-down list next to **Set # Selected To:** on the lower portion of the screen. Click updated **Set # Selected To:**.

Click **Apply** on the very bottom of the page. To undo changes (before clicking **Apply**) click **Reset Changes**.

3.5. Users Tab

To add a new user, click **Add User**. Enter login name and password, and re-enter password to confirm in the pop-up box. Click **OK**.




To inactivate users, simply uncheck the box next to the users' names. Check *Hide Inactive* box to take inactive users out of the users list. This option does not delete users' names from the list, only hides them.

Users tab has four sub-tabs: Basic, Permissions, Groups, and User Plants.

3.5.1. Basic Sub-Tab

To change login password, go to the *Basic* tab and click **Change Password**. Enter password. Reenter password in the next box to confirm. Click **OK**.




3.5.2. Permissions Sub-Tab

To set or change permissions for a user, click on *Permissions* sub-tab.


Click on the user's name to highlight it.

To change permission settings individually for each category, click on each row in the *Grant* column, a drop-down list of permission settings will occur. Make the desired selection.

Click **Reset Changes** to undo changes. Click **Apply** to save changes.




To change permissions for multiple categories at once, hold down the Ctrl key and select desired rows.






The first button directly below the permissions list will update the number of rows selected. Click on the second button, a drop-down list with permission settings will appear. Make the desired selection, then click the updated **Set # Selected To:** button to left of it.



Click **Reset Changes** to undo changes. Click **Apply** to save changes.


3.5.3. Groups Sub-Tab

To assign users to groups, simply click on the user name, then select the group name(s). Click **Apply**.



3.5.4. User Plants Sub-Tab

To assign users to plants, simply click on the user name, then select the plant name(s). Click **Apply**.



3.6. Printers Tab

Report Type drop-down list includes reports with various dimensions used by the plant. Select one.

Printer Name is the printer name (identification code) used to print the report. Select the appropriate printer used for the report dimensions selected.

The *Driver* and *Port* information will be filled automatically.

Click **Apply**. To undo any changes made before clicking **Apply**, use the **Reset Changes** button.


If there are multiple report types, follow the steps above to set appropriate Printer Name to other reports. Click **Apply** after each setting.

3.7. Settings Tab

This screen will be used initially to customize the software to the needs of the company. If need be, simply make changes by typing them into the *Value* column.


3.8. Locations Tab

To add locations or make changes to the location hierarchy, first select the appropriate plant from the *Plant* drop-down list.



Highlight an existing location from the *Location Hierarchy*, fill out the *Location Information* for the new location. Fill out the *Loc. Code* and *Description*, check the box if the new sight is testable, and select the appropriate product proximity and process zone from the drop-down lists. The *Severity Class* is filled automatically based on the process zone selected. Click on the arrow next to **Insert** and specify where in the hierarchy the new location belongs: parent or child of the location highlighted.

Note: Severity classes are assigned to product proximity and process zone combinations in the *Analysis Management: Severity Classes* screen.



The image shows a software interface with two main panels. The top panel, titled "Location Information", contains the following fields and controls:

- Loc. Code**: A text input field.
- (Optional) Description**: A text input field.
- Testable**: A checkbox that is currently checked.
- Product Proximity**: A dropdown menu with a pencil icon to its right.
- Process Zone**: A dropdown menu with a pencil icon to its right.
- Severity Class**: A text input field containing the word "High".
- Create new location...**: A button with a green plus icon and a dropdown arrow.
- Selected Location...**: A section containing four buttons: "Load" (with a green double arrow icon), "Update" (with a blue floppy disk icon), "Cut" (with a blue scissors icon), and "Copy" (with a blue document icon).
- Below the "Selected Location..." section are two more buttons: "Paste" (with a dropdown arrow) and "Delete" (with a red X icon).

The bottom panel, titled "Global Configuration", contains:

- Delimiter character(s)**: A text input field containing a hyphen "-".

A small icon of a computer monitor is visible in the bottom left corner of the interface.

Use **Load** button to bring up the *Location Information* of the location selected from the *Location Hierarchy*.

Use **Cut**, **Copy**, and **Paste** buttons to adjust the *Location Hierarchy*. Select a location, click **Cut**. Click on the arrow next to the **Paste** button, specify if the cut location needs to be pasted below or above the location selected in the *Location Hierarchy*.

To delete a location, select a location from the *Location Hierarchy* then click **Delete**.


Click on **Update** to see the changes under the *Location Hierarchy*. This option does not save the changes made.

Click **Apply** to save changes. To undo a change (before clicking *Apply*) use the **Reset Changes** button.

The pencil buttons next to the *Product Proximity* and *Process Zone* are shortcuts to the *Severity Classes* tab under the *Analysis Management* window. Changes to the product proximity, process zones, and severity class combinations can be made in this screen.

3.9. Key Shortcuts Tab

This page is for reference and to create, change, and delete shortcuts. For example, setting F5 key as a shortcut for opening the *Analysis Management* screen enables the user to press the F5 key (while the Presage Software is open) to open the screen instead of using the mouse to click on the *Analysis Management* under the *Main Menu*. Any key can be used as a shortcut.



To set a shortcut, select one of the options under the column *Name* and click **Set**. A new pop-up window comes up. Press the key to be used as a shortcut, then click **OK**.


To delete a shortcut, select the one of the options under the column *Name* and click **Clear Key**.

To use default shortcuts, select one of the options under *Name* and click **Use Defaults**.

For any changes made, click **Apply** to save. To undo changes (before clicking **Apply**), use the **Reset Changes** button.


4. Data Management Screen

This screen is used to import data into Presage.




5. Investigation Screen

To send an analysis for investigation, first find the failed analysis from the *Work Order* screen or *Sampling History*; right-click on it, and select **Begin investigation from # selected samples**.






User is automatically taken to the *Investigation* Screen.



Right-click on the analysis on this page to see multiple action options: view the location visualization or send to the work order. There is also a **Send to WO** button on the lower, right-hand portion of the screen. Click on the arrow next to the **Send to WO** button, choose between *Send to Existing WO* and *Send to New WO*.



Drag the bar to the right to see the column titled *Comments* and *Findings*. Type in the findings and corrective actions taken in the space provided.

If there are more fails related to this investigation, send the next failed analysis to the investigation screen. This fail will be listed in the *Investigation Detail* (lower) box. Make sure to send the initial analysis listed in the upper box in the *Investigation Screen* to the WO for a retest.

5.1. Configure Investigations

5.1.1. Constraints Sub-Tab


This screen is currently under construction.

5.1.2. Alerts and Notifications Sub-Tab

Use this screen to set up e-mail alerts which can be sent from the 'Work Order' screen by checking the box next to *Send Notification* in the 'Work Order' screen and clicking on *Save*.

First, select a plant. Second, enter the e-mail addresses for the recipients. Third, select which results should be included in the e-mail: fails, marginals, and/or acceptable results.


The template can be customized to each company's needs, simply contact a Presage representative with requests.



6. Locations Screen

This page is for creating a visual map of the plant or location.

First, select the plant from the *Plant* drop-down list in the upper right-hand portion of the screen.




To upload a blueprint of the plant, click on **Load Picture**, locate and select the file, then click **Open**.

Click **Zoom to Fit** to make the entire map visible.

Uncheck the *Show gridlines* box to get rid of the background gridlines.

To map a location, select a location from the *Location Hierarchy*, then click on the button with a box on the *Mode* pop-up window. Click and drag cursor to place the box where appropriate. Now, the box and the location identification name are connected, and once the location identification name or code is selected, the coinciding box on the map will be highlighted automatically.



To add text, click on the box with the letter A then click and drag on the canvas to see the cursor appear.

Use the magnifying glass button on the *Mode* pop-up window to zoom in the picture or diagram.


Use the arrow button to rearrange and re-size boxes. Use the hand button to drag the canvas.

Tip: When several small squares are laid out very close to each other, it is easier to select one particular box by clicking on any box within the group and pressing Tab on keyboard to select other boxes within that area.

For guidelines on *Location Information*, refer to *Configuration: Locations*.

7. Product Management Screen

Product Management is where ingredients and products can be managed. In this screen, users can keep an up-to-date global list of products and assign and remove products to and from each plant.



8. Report Viewer

The report viewer previews reports before they are printed, making it easy to find the correct report and to make necessary adjustments.

8.1. Report Selector

The *Report Name* lists the names of available reports for each report type. Click on a report. *Description* area gets automatically filled out with information about the report.

The *Parameters* section shows what extra information Presage needs to generate a report. Some reports require date ranges, while others require a user to type in information.

Click on the appropriate parameter, then input values into the box provided. For date parameter, click on the drop-down menu then the calendar icon to select the necessary date.

Click on the *Preview* button to open the *Crystal Reports Viewer* tab which displays the preview of the report before printing.

Click *Print* to print the selected report.

8.2. Crystal Reports Viewer

The *Export Report* button exports the report into Excel format. On the pop-up window, select *Microsoft Excel 97-2000(XLS)* from the *Format* drop-down menu. Make sure the *Destination* drop-down menu has *Disk file* selected. Click *OK*. Another window will pop-up. In the new window, define the formatting options for the Excel file. Click *OK* and specify where to save the report.







Reports can also be exported in other formats such as CSV (comma separated value) and XML (eXtensible Markup Language).

The *Print Report* button prints the report from the *Crystal Reports Viewer*.

Toggle Group Tree button will slide over the report preview and open up a new section of the screen. The new section will show how many items are in each group and subgroup of data (if it is applicable for the report).

The *Crystal Reports Viewer* has a group of navigation buttons near the top of the page.

Image	Button
	Go to the first page
	Go to the previous page

 	Go to the next page
 	Go to the last page
 	Go to page x of y

9. Sampling Screen


This screen is used as a search engine which makes it easy to find past sampling information.

Check the box(es) next to the analysis name(s) which pops-up the list of options: check the appropriate boxes for options; plug in the variables in the boxes provided: date range, plants, location, user, proximity and/or zone.

Click **Search**.

Once historical data is pulled up on the sampling screen, four things can be done:


1) Sent to the scatter graph for trending as explained above. 2) Sent to the Location Visualization which shows the sample acceptability over the plant layout. 3) Sent to the Graphing to do more complex charting. 4) Exported to a .csv file.



Click on the column title on the table, to order the data by a certain criteria. For example, click on the *Performed On* to see results in chronological order of date performed.

To Trend Data:

Click the *Graphing* tab for a scatter graph for the criteria selected.





Select and change criteria for the graph using the same panel as used with sampling history on the upper portion of the screen.

Hover over points to view details such as date, time, analysis, analysis option, location, etc.

Hover over a label in the map key to highlight the group of points which belong to that label.


To share the graph with a third party, click on the copy icon next to the browser link box on the mid-right side of the screen (it is the second button next to the browser link box). Then paste the link.

Check the **Show Thresholds** button to see the threshold lines on the graph.

9.1. Visualization

To see Location Visualization:

First, enter criteria and pull up historical data. Right-click on one of the rows. A dialog with multiple options pops-up. Click **send all # results to the locations screen for visualization**.



On the left display, locations are highlighted in green and red. Green locations are ones with acceptable counts, and red locations are ones with fail counts.

Highlight a row on the right and the corresponding location will be highlighted automatically on the left display and vice versa.

Order the *Date Sampled* column by clicking on the column title **When**. Click the first row and press the down-arrow on the keyboard to see a play-through of where the samples were taken in a chronological order.

Note: Every row can be ordered by clicking the column title.

Check the box **Color locations based on selected samples only**, and select a group of rows to see the colors only for those selected.




9.2. Graphing


To do more complex graphing:

Once historical data is pulled up on the sampling screen, this data can be graphed in multiple ways, in addition to the scatter graph.


First, enter criteria and pull up historical data. Right-click on one of the rows. A dialogue with multiple options pops-up. Click **export all rows to graph**.



The following dialog pops-up. Make selections, then click **OK**. The graph will open in a new browser page.



These graphs are saved in the local computer and cannot be shared with third-parties.



9.3. Export Data

Please, contact a Presage representative for exports.

10. Work Orders


10.1. Edit Work Order Tab

Use this page to edit work orders, add additional analyses to work orders, and enter data collected for analyses.

Click on the work order number in the *Home* screen to pull it up in this page. There is also a way to search for work orders in this page by clicking on the blank box under *Work Order:*. On the pop-up window, enter the work order number. Click **OK**. The work order will show up under the *Sample Detail* tab.

Now, enter the results for the analyses under the columns *Option Value*. There may be more than one depending on how many analysis options were assigned under *Analysis Management*. Make sure the screen is scrolled all the way to the right to see all analyses options and/or questions. Grayed out boxes are blocked and there is no option to enter data there.

To collapse some columns, uncheck the appropriate boxes on the lower portion of *Sample Detail* tab. The *Show analysis and option values* box collapses columns titled *Option Name* and *Option Value*. The *Show collection detail* box collapses columns titled *Performed On* and *Performed By*. The *Show testing detail* box collapses columns titled *Testing Began*, *Testing Ended* and *Tested By*.



Under *Sample Detail*, click on the box under *Status* to change the status of the work order. Select the appropriate status title from the drop-down menu: *Opened*, *Sampled*, *Closed*, and *Cancelled*. Set the appropriate date and time the analysis was performed in the next column. Then select the appropriate user who performed the test in the column *Performed By*. *Tag #* is assigned automatically to analysis, but there is an option to change it; simply click on the box to make the change. Make sure to select a unique number for each analysis. For those who are using the scanner, these columns are automatically updated when the scan is performed.

To add an analysis, click **Add**, a new highlighted row will pop up. Start entering data for the new analysis. To remove an analysis, select the row (use CTRL key to select multiple rows), then click **Remove** on the lower portion of the screen.

The **Set Status...** button on the lower left corner of the screen allows users to change the status of groups of analyses simultaneously. Select **All Items...** or **Only Selected Items...**, then **To Cancelled** to invalidate the analyses in the Work Order. Invalidated samples do not show under *Sample History*.


To assign analyses to groups, select the analysis name from the list, then select the group name from the drop-down menu under *Assign Sampling to Group*:. Click **Distribute**.

To save all changes, click **Save**. There are options to print and close work order once changes are saved; if applicable, click on the check-box.

When printing reports or tags, check the **&Print** box when saving the work order. A window will pop-up with list of reports, choose the appropriate report. On the next window, choose the appropriate printer and click **OK**.

10.2.Sample Schedule

This screen lists all the analysis tests which need to be performed. First column lists the analyses which are created in the *Analysis Management* screen. Second column lists the plant that particular analysis is scheduled to be tested in. Next column is the site the analysis is assigned to be sampled from. *Schedule Name* is the interval period the test is performed in; it is created in the *Analysis Management*. *Date scheduled* is the date the analysis was created and *Time Remaining* is the number of days left to perform the test.



Choose the plant from the drop-down menu on the right.

Choose the date range to view analyses scheduled for specific dates. For custom date range, select *Custom* from the drop-down menu. Then set the *From:* and *To:* dates from the calendar icon or type it in. Otherwise, leave the date range as is, as this schedule is ordered by the least number of days remaining to be performed.

Check the **Show samples already on Wos** box to see all analyses, ones already assigned to a work order in addition to ones not assigned yet. The work order numbers will show up in the second to last column, and work order status will show up in the last column.

The screenshot shows a software window titled 'Edit Work Order' with a 'Sample Schedule' tab. At the top, there are controls for 'Date Range' (set to 'Custom'), 'From' and 'To' dates (both set to 5/23/2013), and checkboxes for 'Show samples already on WOs' and 'Show samples past due'. A 'Plant' dropdown is set to 'LNK - LINCOLN PLANT'. Below these controls is a table with the following data:

Analysis	Plant	Location	Schedule Name	Date Scheduled	Time Scheduled	Time Remaining
PreOp	LINCOLN PLANT	ST1-2	Semi-Monthly P...	5/24/2013	9:00 AM	21 hours, 1 minute

A red rectangular box highlights the 'Time Remaining' column header and its value in the data row.

Click on the column titles to put the list in order based on the selected category. For example, click on the *Analysis* column title to order the list alphabetically by analysis names; click on *Location* to order the list alphabetically based on the name of locations, etc.


This list is color-coded. These colors are the default and can be customized.

Blue samples need to be completed within 24 hours.

Green samples have more than a day left to be completed.

Grey samples are past due.

Note: Past due samples will not automatically move to current schedule. They will have to be manually selected and sent to work orders. Set the date range to the past to see past due orders.



To assign random sights to a work order, click on **Choose Random**. A new window will pop-up. Specify the number of locations the system should pick out randomly under the *Choose* column.


Click **Choose Items**. Random analyses will be selected in the Sample Schedule.

Now, click **Send to WO**. The highlighted analyses from *Sample Schedule* will now show up under the *Edit Work Order* tab.

Make sure the correct *Date Scheduled* and *Plant* are selected.

Click **Save**.

Now, the work order is created and is listed under *Work Orders* in the *Home Screen*.



There is also a way to make manual selections. Press on the *Ctrl* key while making multiple selections on the schedule. Click **Send to WO**, then **Save**. Click **Clear** (before **Save**) to put samples back to the schedule.

If multiple analyses are assigned to one location, the system can grab all analyses tied to that location to the work order. So, if a user specifies ten locations to be randomly picked out from the *Sample Schedule*, there may be more than ten analyses in that work order depending on how many analyses are assigned to each location chosen.

To achieve that, simply make selections from the schedule, click on the arrow-down button next to **Send to WO**, check the box next to the pop-up message **Perform all scheduled analyses for each chosen location**, then click **Send to WO**.

11. How to insert the iPod in the Linea Pro scanner.



11.1. Install the ITrack Mobile App on the iPod.

Click on the **App Store** icon on the iPod.

Type in *isoftdata*, and swipe from right to left until the following app shows up:





**The button may look different.*

Click on the cloud button to install the app.

The following icon pops-up. Click **Use Existing Apple ID**.



Enter the Apple ID and password. (A Presage representative will provide you with this information. It will be used across all plants in your companies, unless requested otherwise.)




In the Home Screen, find the ITrackLX Mobile app icon and click on it.





Enter the Server URL provided by the Presage representative.



Login with Presage username and password and click **Log in!**.

11.2. Perform Sample

To time-stamp when a sample was taken, click on the first button **Perform Sample**.

Press the button on the right side of the Linea Pro.

Place the red ray on the barcode. Hold until a sound comes out validating that the barcode was read.

A notification will appear confirming that the sample was marked.

Click **Add Image** to attach an image to the sample.

Click **Done** to go back to the main menu.

The above steps will automatically update the work order in Presage.

- 1) The *Status* will change from *Open* to *Sampled*.
- 2) *Performed On* will be filled with the scan date and time.
- 3) *Performed By* will be filled with the user who logged in to the mobile app.

11.3. Incubate Sample

To time-stamp when a test is going in and out of the incubator, click on the second button **Incubate Sample**.

Press the button on the right side of the Linea Pro.

Place the red ray on the barcode. Hold until a sound comes out validating that the barcode was read.

A notification will appear confirming that the sample was marked.

Click **Add Image** to attach an image to the sample.

Click **Done** to go back to the main menu.

The above steps will automatically update the work order in Presage.

- 1) *Testing Began* will be filled with the date and time the barcode was scanned the first time.
- 2) *Testing Ended* will be filled with the date and time the barcode was scanned the second time.
- 3) *Tested By* will be filled with the user who logged in to the mobile app.

11.4. Record Analysis

To enter the result for a sample, click on the third button **Record Analysis**.

Press the button on the right side of the Linea Pro.

Place the red ray on the barcode. Hold until it beeps and sample appears on the screen.

All the options for the sample will show up. Select or type in the results.

Click **Done** to save and go back to the main menu. Click **Cancel** to delete changes and go back to the main menu.

The above steps will automatically update the work order in Presage.

12. How to:

How to create a work order.

1. Click on the **Work Order** screen under *Main Menu*.
2. Click the **Add** button on the lower-left of the *Edit Work Order* tab.
3. Once the row is selected, click on the cell below *Analysis* to choose an analysis name.
4. Click on the cell below *Location* to choose a location name.

Note: If there are multiple analyses to be entered, click the **Add** button multiple times to make the rows editable.

Note: Make sure to remove any blank editable rows before saving the work order. Select the row to be deleted, and click **Remove**.


5. Check the appropriate boxes next to printing options.

Remember: **&Print** prints letter size reports.

&Print Tags has two options: with Barcode and without Barcode. Usually it's a 1.1 x 2.4 label, but it can be customized to your needs.

&Print Test Tags prints very thin labels with barcodes. These are meant for sticking on petri dishes and petri films.

6. Click **Save**.



How to hide unnecessary columns from the work order screen.

Uncheck the appropriate boxes on the lower portion of work order screen.

The **Show analysis and option values** box collapses option name columns.

The **Show collection detail** box collapses columns titled Performed On and Performed By.

The **Show testing detail** box collapses columns titled Testing Began, Testing Ended and Tested By.

How to open a closed work order.

- 1) Pull up the work order in the work order screen.
- 2) Click on the **Set Status** button on the lower left, select **Re-open closed work order**.

How to edit a closed sample.

- 1) Pull up the work order in the work order screen.
- 2) Click on the **Set Status** button on the lower left, select **Only Selected Items...** or **All items...**, then select **To Sampled**.
- 3) Click **Re-Open WO** on the pop-up dialog.


How to invalidate a sample.

In the work order screen, select the row, click on the arrow next to **Set Status...** on the bottom-left corner. Hover over **Only Selected Items** or **All Items** then click on **To Cancelled....**


Invalidated samples do not show in the Sample History.

How to send an analysis for investigation.


To send an analysis for investigation, select the sample analysis from the *Work Order* screen or *Sampling History*; right-click on it, and select **Begin investigation from # selected samples**.



User is automatically taken to the *Investigation* Screen.



Right-click on the analysis on this page to see multiple action options: view the location visualization or send to the work order. There is also a **Send to WO** button on the lower, right-hand portion of the screen. Click on the arrow next to the **Send to WO** button, choose between *Send to Existing WO* and *Send to New WO*.



Drag the bar to the right to see the column titled *Comments* and *Findings*. Type in the findings and corrective actions taken in the space provided.

If there are more fails related to this investigation, send the next failed analysis to the investigation screen. This fail will be listed in the *Investigation Detail* (lower) box. Make sure to send the initial analysis listed in the upper box in the *Investigation Screen* to the WO for each retest.


How to set up a sample.

- 1) Click on **Analysis Management** under the *Main menu*. Stay on the first tab *Analyses*.
- 2) To create a new analysis, click on the box *<Enter Analysis Name>*. Enter the analysis name on the box provided. In the next column, type in the category name. This column is used to group analyses. *Testing Time* is for lab use. Type in the number hours the sample must rest in the incubator.
- 3) Under *Option*, enter what is being tested for; for example: plate count, sample area, detect-ability, etc.

When appropriate, enter the unit in the next column; for example: cfu, ppm, etc.

The **Informational** box is for options which don't need to be graphed. Check if appropriate.

Then choose the *Value Type* from the drop-down list provided. For example, plate count is measured in number, and qualitative analyses are measured in boolean: positive or negative. To create a drop-down list of various options to select from, choose *Choice*.



An analysis disappeared from Analyses screen, and I did not delete it.

If an analysis has a global setting, then it is used by all plants within a company. If one plant unchecks the global setting, then the analysis is hidden for all other plants except for the plant that set the analysis to be plant-specific.

Note: This does not mean the analysis and all data attached to it is deleted for all other plants. It is simply hidden.

Make sure another plant has not unchecked the global setting in the *Analysis Management-Analyses* screen.

If not, call the customer support rep at 1.800.309.1704.


How to set up a threshold.

- 1) Click on **Analysis Management** under the *Main Menu*. Stay on the first tab *Analyses*.
- 2) Select an analysis and a specific option you want to set thresholds or specifications to. Set thresholds for each option separately.
- 3) Select the severity from the drop-down list. To set up severity classes, refer to *Analysis Management: Severity Classes*.
- 4) For *Threshold*, select an option from the drop-down list. Marginal values are those which are acceptable, but are nearing the unacceptable boundary. Unacceptable values represent failure. Invalid values signify an error.

5) *Constraint Type* sets the minimum and maximum boundary value.

6) Choose the appropriate *Choice/Boundary Value*, or fill it in manually. For example, to set 10 cfu as a maximum acceptable limit for an analysis, choose *Unacceptable when* under *Thresholds*; *Above* under *Constraint Type*; and type in *10* under *Choice/Boundary*.

7) Click **Apply** to save changes.



How to make thresholds plant-specific for global analyses.

Select a plant-specific *Severity* with each threshold. **Default** *Severities* are global and will show for all plants.

Analysis Management

Analyses | Manage Schedules | Analysis Scheduling | Severity Classes

Plant: LNK - LINCOLN PLANT

Name	Global	Category	Group...	Testing Ti...
TPC	<input checked="" type="checkbox"/>	Product Testing	<input checked="" type="checkbox"/>	48
Enterobacteriaceae	<input checked="" type="checkbox"/>	Environmental	<input checked="" type="checkbox"/>	24
e-coli	<input checked="" type="checkbox"/>	Product Testing	<input checked="" type="checkbox"/>	24
Consumer Complaints	<input checked="" type="checkbox"/>	Complaints	<input checked="" type="checkbox"/>	
Trap Inspection	<input checked="" type="checkbox"/>	Maintenance	<input checked="" type="checkbox"/>	
Wash Concentration	<input checked="" type="checkbox"/>	Maintenance	<input checked="" type="checkbox"/>	
Allergen	<input checked="" type="checkbox"/>	Product Testing	<input checked="" type="checkbox"/>	
Product Mapping	<input checked="" type="checkbox"/>	Environmental	<input checked="" type="checkbox"/>	
Shelf Life	<input checked="" type="checkbox"/>	Product Testing	<input checked="" type="checkbox"/>	
ATP	<input checked="" type="checkbox"/>	Environmental	<input checked="" type="checkbox"/>	
QC general	<input checked="" type="checkbox"/>	QA	<input checked="" type="checkbox"/>	
Internal Audits	<input checked="" type="checkbox"/>	Audits	<input checked="" type="checkbox"/>	
Finished Product Yes...	<input checked="" type="checkbox"/>	Product Testing	<input checked="" type="checkbox"/>	
<Enter Analysis Na...	<input checked="" type="checkbox"/>		<input checked="" type="checkbox"/>	

Analysis options/questions

Option	Unit	Default Value	Value Type	Required	Ir
ProdDatetime			Datetime	<input type="checkbox"/>	
Coliform			Number	<input type="checkbox"/>	
e-coli			Choice	<input type="checkbox"/>	
pH			Number	<input type="checkbox"/>	
HPC			Number	<input type="checkbox"/>	
Yeast/Mold			Number	<input type="checkbox"/>	
<Enter Analysis ...			Number	<input type="checkbox"/>	

Option choices/boundaries

Severity	Threshold	Constraint T...	Choice/Boundary
High	Unacceptable when	Above	10
Default	Marginal when	Above	<Enter Option...

How to set up a schedule.

- 1) First select the plant from the drop-down list.
- 2) Under the column title *Analysis*, select the analysis type from the drop-down list. This list includes all analyses created under *Analysis Management: Analyses*.
- 3) Select the schedule period interval from the drop-down list under *Schedule*. This list is created under *Analysis Management: Manage Schedules*.
- 4) The *Location Count* is automatically filled out when the locations are selected in the next step.

Analysis/Schedule:

Analysis	Schedule	Loc Count	
TPC	Semi-Monthly Pre...	103	<input type="checkbox"/>
Enterobacteriaceae	Semi-Monthly Pre...	1	<input type="checkbox"/>
Sanitation 1	Sanitation Pre-Op	1	<input type="checkbox"/>
		0	<input type="checkbox"/>

- 5) The right side of the screen is for choosing locations to perform the analyses on. Select one analysis under *Analysis/Schedule*. Then check the box next to the locations which need to be sampled.
- 6) For a quick look-up, type in the first letter or number of the display code for the location. The corresponding locations will show up in the display below. The list can also be filtered by severity class. Choose the severity class from the drop-down; the list in the display will narrow down to locations with the severity class specified.

Add/Remove locations quick lookup:

Severity Class Filter

7) Use the **Check...** button to select groups of locations. Use the **Uncheck...** button to deselect groups of locations. Click **Add** (on the middle of the screen, not the one on the lower portion of the screen) to add the selected locations to the list of locations to be sampled. Click **Remove** to take out locations from the list.

Code	Location	Description
<input type="checkbox"/> ST1	ST1	Stage 1
<input checked="" type="checkbox"/> 1	ST1-1	Stage 1-1
<input checked="" type="checkbox"/> 2	ST1-2	Stage 1-2
<input checked="" type="checkbox"/> 3	ST1-3	Stage 1-3
<input checked="" type="checkbox"/> 4	ST1-4	Stage 1-4
<input type="checkbox"/> ST2	ST2	Stage 2
<input checked="" type="checkbox"/> 5	ST2-5	Stage 2-5
<input checked="" type="checkbox"/> 6	ST2-6	Stage 2-6
<input type="checkbox"/> ST3	ST3	Stage 3

Locations to be sampled:

Location	Plant
ST1-2	LNK - LINCOLN PLANT
ST1-3	LNK - LINCOLN PLANT
ST1-4	LNK - LINCOLN PLANT
ST8-94	LNK - LINCOLN PLANT
ST8-95	LNK - LINCOLN PLANT
ST8-96	LNK - LINCOLN PLANT
ST8-97	LNK - LINCOLN PLANT
ST8-98	LNK - LINCOLN PLANT
ST8-99	LNK - LINCOLN PLANT


8) To add locations one at a time, use the **Add** button on the lower portion of the screen under *Locations to be sampled*. To remove a location, select a location from the list and click **Remove**. To remove multiple locations, hold down the **Ctrl** button on keyboard while selecting multiple locations, then click **Remove**.

9) Click **Apply** to save any changes and use the **Refresh/Reset** button to remove any unsaved changes.

How to set up a severity/proximity/zone matrix.

1) Click **Analysis Management** under the *Main Menu*.

2) Click **Severity Classes**.



3) Give a unique names to all severity classes. Severity classes are based on limits which are set based on product proximities and process zones. Examples are High, Medium, and Low; or 10, 100, 500. Use terminology or numbers which best fit your system.

4) Click on the box under *Default* (may have to expand the column to see the full column title), to set a default. Default severity class is usually the one used the most often.

5) *Product proximity* refers to contact or non-contact surfaces. To add new product proximities, click on <New Proximity>. Type in the new proximity title. To remove product proximity, click on the red box with a white X to the right of each line. To rename, select the appropriate row and type in the new proximity name.

6) *Process zones* are different operational zones such as Chilling or Packaging. To add new zones, click on <New Zone> (It may be necessary to scroll down to see this row.) To delete a zone, click on the red box with white X to the right of each line. To rename, select the appropriate row and type in the new zone name.

7) To assign severity classes to zone/proximity combinations, click on one severity class under *Severity Classes* on the left display to highlight it, then check the boxes where the selected severity class is appropriate on the right display.

Process Zone	Contact	Non-Contact	Environmental
Zone 1	<input checked="" type="checkbox"/> High	<input type="checkbox"/> Default	<input type="checkbox"/> Low
Zone 2	<input type="checkbox"/> Medium	<input type="checkbox"/> Default	<input type="checkbox"/> Default
Zone 3	<input checked="" type="checkbox"/> High	<input type="checkbox"/> Default	<input type="checkbox"/> Low
Zone 4	<input type="checkbox"/> Default	<input type="checkbox"/> Default	<input type="checkbox"/> Default
Zone 5	<input checked="" type="checkbox"/> High	<input type="checkbox"/> Medium	<input type="checkbox"/> Default
Zone 6	<input type="checkbox"/> Default	<input type="checkbox"/> Default	<input type="checkbox"/> Low
Zone 7	<input checked="" type="checkbox"/> High	<input type="checkbox"/> Medium	<input type="checkbox"/> Default
Zone 8	<input type="checkbox"/> Medium	<input type="checkbox"/> Default	<input type="checkbox"/> Low
Zone 9	<input checked="" type="checkbox"/> High	<input type="checkbox"/> Medium	<input type="checkbox"/> Default

8) Now, click the next severity class to highlight it, then check the boxes where this severity class is appropriate on the right display. Not all boxes have to be checked.

9) Click **Apply** to save the changes.

How to add a plant.

- 1) Click **Configuration** under the *Main Menu*.
- 2) Click **Plants**.
- 3) Click **New Plant**.
- 4) For *Plan #*, select *New* from the drop-down list.
- 5) If the new plant is private, check the box provided. With this option, users must log in to have access to the new plant.
- 6) The *Display Code* will appear in the *Plant #* drop-down list after the new plant has been created. Choose a unique code which will best identify the plant.
- 7) *Company Name* refers to the business name a particular plant is associated with.
- 6) Fill in all information in the spaces provided.


How to delete a plant.

- 1) Click **Configuration** under the *Main Menu*.
- 2) Click **Plants**.
- 3) To delete a plant, choose the plant from the *Plant #* drop-down list.
- 4) Click **Delete Plant**.

How to add a user.

- 1) Click **Configuration** under the *Main Menu*.
- 2) Click **Users**.
- 3) Click **Add User**.
- 4) Enter login name and password, and re-enter password to confirm in the pop-up box.
- 5) Click **OK**.

Note: Password is case sensitive, but the username is NOT case sensitive. Choose a password that is easy to remember but difficult for others to guess. A good password has a combination of both letters and numbers.



How to delete a user.

In Presage, there is no way to completely delete a user, however users can be inactivated and hidden from the users list.

- 1) Click **Configuration** under the *Main Menu*.
- 2) Click **Users**.
- 3) Uncheck the box next to the users' names.
- 4) Check *Hide Inactive* box to take inactive users out of the users list.

How to change a user password.

- Note: Password is case sensitive. Choose a password that is easy to remember but difficult for others to guess. A good password has a combination of both letters and numbers.

- 1) Click on **Configuration** under the *Main Menu*.
- 2) click on the **Groups** tab.
- 3) Click **Add Group**.
- 4) Type in the new group name.
- 5) Click **OK**, then **Apply**.
- 6) Check the box under *Show on WOs*, if this group will perform tasks in the work orders.

User Count shows the number of users in the group.



How to set permissions to groups.

1) Click on **Configuration** under the *Main Menu*.

2) Click on the **Groups** tab.

Permission settings are *none* (not permissible), *plant* (permissible only for this plant), and *global* (permissible for every plant).

3) Scroll to the right to see complete descriptions for each permission.

4) Click on the group name to highlight it.

5a) To change permission levels individually for each category, click on each row in the *Grant* column, a drop-down list of permission levels will appear. Make a selection.

5b) To change permissions for multiple categories at once, hold down the **Ctrl** key on the keyboard and select desired rows. The first button directly below the permissions list will update the number of rows selected. Click on the blank button next to it, a drop-down list with permission levels will appear. Make a selection, then click the updated **Set # Selected To:** button to left of it.

6) Click **Apply** to save changes.

Note: Users get the best of user and group permission levels. For example, if the user is assigned global permission on the user level, but plant permission for the group level, then that user gets the global level for that specific action.


How to assign users to groups.

1) Click on **Configuration** under the *Main Menu* Then click on the **Users** tab.

2) Click on **Groups** which is the third sub-tab under the **Users** tab.

3) Click on the user's name to highlight it, then click on the group name.

4) Click **Apply**.



How to set permission for each user.

Go to **Configuration**, then **Users**.


Click on the user's name to highlight it, then click on *Permissions* sub-tab.

To change permission settings individually for each category, click on each row in the *Grant* column, a drop-down list of permission settings will appear. Make the desired selection.


Permission settings are *none* (not permissible), *plant* (permissible only for this plant), and *global* (permissible for every plant).

Scroll to the right to see complete descriptions for each permission.


Click **Apply** to save changes.



To change permissions for multiple categories at once, hold down the **Ctrl** key and select desired rows or **Ctrl-A** to select all.



The first button directly below the permissions list will update the number of rows selected. Click on the second button, a drop-down list with permission settings will appear. Make the desired selection, then click the updated **Set # Selected To:** button to left of it.




Click **Apply** to save changes.

User's user permission is different from his/her group permission.

In this case, the user is given the best of both permission levels. For example, if the user permission for deleting data is set to plant-level and the group (which the user is assigned to) permission for the same task is set to global-level, then the user gets the global permission.

How to assign users to plants.

- 1) Click on **Configuration** under the *Main Menu*.
- 2) Click **Users** tab then **User Plants** sub-tab.
- 3) Click on the user name, then select the plant name(s).
- 4) Click **Apply**.





How to configure printer settings.

1) Click on **Configuration** under the *Main Menu*.

2) Click on the **Printers** tab.

3) *Report Type* drop-down list includes reports with various dimensions used by the plant. Select one.

4) *Printer Name* is the printer name (identification code) used to print the report. Select the appropriate printer used for the report dimensions selected.

The Driver and Port information will be filled automatically.

5) Click **Apply**.

If there are multiple report types, follow the steps above to attach a report type to a printer name. Click Apply after each setting.

How to print labels.


After adding all the samples in the work order, check the appropriate **&Print** box. Then click **Save**.

&Print prints letter size reports. **&Print Tags** has two options: with Barcode and without Barcode. Usually it's a 1.1 x 2.4 label, but it can be customized to your needs. **&Print Test Tags** prints very thin labels with barcodes. These are meant for sticking on petri dishes and petri films.

The screenshot shows a software interface with a table of data and a control panel. The table has two columns with numerical values: 304.3, 302.8, 303.5, and 300.3 in the first column, and 0.08, 0.15, 0.15, and 0.13 in the second column. Below the table is a control panel with an 'Attachments' button, a 'Clear' button, a 'Save' button with a dropdown arrow, and two checkboxes: 'Close WO' and 'Send Notification'. To the right of these is a red-bordered box containing three checkboxes: '& Print', '& Print Tags', and '& Print Test Tags'. At the bottom of the interface, there is a status bar with a green checkmark icon, the text 'Updater: (2:26 PM) software up to date', 'Current User: support @ 7', and 'NUM'.

How to print a report.


1. Click **Report Viewer** under *Main Menu*.
2. Select the report type.
3. Select the report name.
4. Enter a parameter value next to the parameter name listed.
5. Click **Preview** to view the report, or click **Print**.



How to narrow down history search using option filters.


1. Click **Sampling** under *Main Menu*.
2. Check the box next to the plant you want to pull up data on.
3. Check the box next to the analysis you want to search for.
4. Check the box next to the option you want to search for.
5. Click on the field under *Result Status* and select a status from the list to narrow down search to a specific result status.

6. Click on the field under *Filter* to type in or select a value to narrow down search. If the option is set up to be a choice, this field will be a choice. If the option is set up to be a date, then select an appropriate date. If the option is set up to be a number, type in a number to filter.
7. Make appropriate Date Range, Location, Proximity, Zone, and User selections.
8. Click **Search**.
9. Data will appear under the *Sample Detail* tab.



How to add/edit/delete locations.

To add locations or make changes to the location hierarchy, first select the appropriate plant from the *Plant* drop-down list.



Highlight an existing location from the *Location Hierarchy*, fill out the *Location Information* for the new location. Fill out the *Loc. Code* and *Description*, check the box if the new sight is testable, and select the appropriate product proximity and process zone from the drop-down lists. The *Severity Class* is filled automatically based on the process zone selected. Click on the arrow next to **Insert** and specify where in the hierarchy the new location belongs: parent or child of the location highlighted.

Note: Severity classes are assigned to product proximity and process zone combinations in the *Analysis Management: Severity Classes* screen.

Location Information

Loc. Code

(Optional) Description

☒ Testable

Product Proximity

Process Zone

Severity Class

High

Create new location...

+ Insert...

Selected Location...

Load Update

Cut Copy

Paste Delete

Global Configuration

Delimiter character(s)

-

Use **Load** button to bring up the *Location Information* of the location selected from the *Location Hierarchy*.

Use **Cut**, **Copy**, and **Paste** buttons to adjust the *Location Hierarchy*. Select a location, click **Cut**. Click on the arrow next to the **Paste** button, specify if the cut location needs to be pasted below or above the location selected in the *Location Hierarchy*.

To delete a location, select a location from the *Location Hierarchy* then click **Delete**.


Click on **Update** to see the changes under the *Location Hierarchy*. This option does not save the changes made.

Click **Apply** to save changes. To undo a change (before clicking *Apply*) use the **Reset Changes** button.

The pencil buttons next to the *Product Proximity* and *Process Zone* are shortcuts to the *Severity Classes* tab under the *Analysis Management* window. Changes to the product proximity, process zones, and severity class combinations can be made in this screen.

How to create a plant map.

First, select the plant from the *Plant* drop-down list in the upper right-hand portion of the screen.




To upload a blueprint of the plant, click on **Load Picture**, locate and select the file, then click **Open**.

Click **Zoom to Fit** to make the entire map visible.

Uncheck the *Show gridlines* box to get rid of the background gridlines.

To map a location, select a location from the *Location Hierarchy*, then click on the button with a box on the *Mode* pop-up window. Click and drag cursor to place the box where appropriate. Now, the box and the location identification name are

connected, and once the location identification name or code is selected, the coinciding box on the map will be highlighted automatically.



To add text, click on the box with the letter A then click and drag on the canvas to see the cursor appear.

Use the magnifying glass button on the *Mode* pop-up window to zoom in the picture or diagram.

Use the arrow button to rearrange and re-size boxes. Use the hand button to drag the canvas.

Tip: When several small squares are laid out very close to each other, it is easier to select one particular box by clicking on any box within the group and pressing Tab on keyboard to select other boxes within that area.

For guidelines on *Location Information*, refer to *Configuration: Locations*.

Why aren't my locations listed in the **Location** drop-down list in the Work Order screen?

When a location is marked as not testable (by unchecking the **Testable** button under *Location Information*), that location will not show up in the Location drop-down list in the work order.

Make sure to check the **Testable** box for all locations that get tested.

For example, in the location below (Loc. Code 002 : Receiving Combo Dumper), the **Testable** box is unchecked.

Locations Visualization

Plant
LNK - LINCOLN PLANT

Location Hierarchy (select an item to add sub-locations)

Code	Description	Testable	Active
Zone 1		No	Yes
001	Conveyor Combo D...	Yes	Yes
002	Receiving Combo D...	No	Yes
003	Guide Rail on West C...	Yes	Yes
004	Guide Rail on E. Con...	Yes	Yes
005	Receiving Product R...	Yes	Yes
006	East Chutes	Yes	Yes
007	West Chutes	Yes	Yes
008	North Hangline	Yes	Yes
009	Incline Conveyor	Yes	Yes
010	South Hang Line	Yes	Yes
011	Transfer Conveyor	Yes	Yes
012	Schackles in Hang Li...	Yes	Yes
013	Schackle Guide	Yes	Yes
014	Coneline	Yes	Yes
015	Wheel	Yes	Yes
016	Table	Yes	Yes
017	Cutting Board	Yes	Yes
018	Product Bin	Yes	Yes
Zone 2		No	Yes
019	Overhead Conveyor	Yes	Yes
020	Overhead Drop Chutes	Yes	Yes
0..	Investigational Site 1	Yes	Yes
021	Overhead Gates	Yes	Yes

Attachments

Location Information

Loc. Code
002

(Optional) Description
Receiving Combo Dumper

☐ Testable ☒ Active

Product Proximity
Contact

Process Zone
Zone 1

Severity Class
High

Create new location...
+ Insert...

Selected Location...
Cut Copy
Paste Delete

☒ Location draw mode (draw box for selected location)
☐ Insertion mode (new box auto creates child location)

Global Configuration
Delimiter character(s)
-

Why aren't all my testable locations not showing in the **Location** drop-down list in the Work Order screen?

The location drop-downs (for speed purposes) only load the first couple hundred items by default. Typing narrows down the list so you can find items that aren't shown.

How to add a product.

1. Click on **Product Management** under the *Main Menu*.
2. Select the appropriate plant from the *Plant:* drop-down menu.
3. Click an existing line in the *Product List*.
4. Click the down-arrow next to the **Add...** button.

5. Select an appropriate position for the new product: above, next to, or below.
6. A new selected line appears in the list. Fill in the appropriate fields: *Name*, *Category*, *Description*, and make sure to check the appropriate boxes: *Product is in use at the plant* and *Active Product*.
7. Click **Apply**.

How to delete a product.

1. Click on **Product Management** under the *Main Menu*.
2. Select the appropriate plant from the *Plant*: drop-down menu.
3. Click an existing line in the *Product List*.
4. Click **Delete**.

Note: Deleting a product means all data attached to it will also be deleted from the historical samples. An alternative is to *inactivate* the product.

5. A pop-up appears warning about the above note. Click **OK**.
6. Click **Apply**.

How to inactivate a product.

1. Click on **Product Management** under the *Main Menu*.
2. Select the appropriate plant from the *Plant*: drop-down menu.
3. Click an existing line in the *Product List*.
4. Uncheck the box next to *Active Product*.
5. Click **Apply**.

Note: To see these inactivated products, check the box next to *Show inactive products* directly below the *Products* tab.

How to assign a product to a plant.

1. Click on **Product Management** under the *Main Menu*.
2. Select the appropriate plant from the *Plant*: drop-down menu.
3. Click an existing line in the *Product List*.
4. Check the box next to *Product is in use at plant "X"*.
5. Click **Apply**.


How to attach an image to a product.

1. Click on **Product Management** under the *Main Menu*.
2. Select the appropriate plant from the *Plant*: drop-down menu.
3. Click an existing line in the *Product List*.
4. Click on the **Attachments** tab.


5. Click **Add**.
6. Select an image and click **Open** on the pop-up.
7. Click **Apply**.

How to apply updates.


When updates have been released, the following dialog pops-up upon opening the Presage application.




Once the following dialog pops-up, click **Apply Updates**.



Then click **Launch ITrack**.



Note: If the updates have not been applied, the notification will appear on the very bottom-right panel of Presage.



Updater: (2:11 PM) No new updates. 1 updates require restart

Double-click on this notification to start the downloading updates process.

My update finished with errors. What do I do?

1. Reset the updates.



2. If the above does not work, call the Presage customer support at 1.800.309.1704.

Note: Presage can be used even if the updates have not gone through by clicking **Close** on the update pop-up dialog and logging in as usual.

How to send samples from the Sample Schedule to the work order.

- 1) Click on **Work Orders** under the *Main Menu*.
- 2) Click on the **Sample Schedule** tab.
- 3) Choose a plant from the drop-down menu on the right.
- 4) Choose the date range to view analyses scheduled for specific dates. For custom date range, select *Custom* from the drop-down menu. Then set the *From:* and *To:* dates from the calendar icon or type it in. Otherwise, leave it to *Today*.
- 5a) Select a sample. To select multiple samples, hold on to **Ctrl** key and select the samples.
- 5b) For random selection, click **Choose Random** on the bottom-right corner. Type in number of samples to pull under the *Choose* column in the appropriate *Schedule* row.
- 6) Click **Send to WO**.

If I add a sample to a work order manually, will it affect the status of the same sample that's in the schedule?

No. If a sample has been added to a work order through the **Add** button instead of the **Send to WO** in the *Sample Schedule* screen, these samples are independent and don't mark the same sample that is scheduled as completed. The scheduled sample will remain in the schedule.

How to view/graph historical data.

- 1) Click on **Sampling** under the *Main Menu*.
- 2) Stay on the *Sample Detail* tab to view data in spreadsheet form. Click on the *Graph* to view data in a scatter graph.
- 3) Select a plant.
- 4) Select an analysis name.
- 5) Check the boxes next to options you want pulled up. Enter filters when needed.
- 6) Plug in the variables in the boxes provided: date range, plants, location, user, proximity and/or zone.
- 7) Click **Search**.

Sampling

Sampling History

Search Criterion

Use	Code	Name
<input checked="" type="checkbox"/>	UNK	LINCOLN PLANT
<input checked="" type="checkbox"/>	CHI	CHICAGO PLANT
<input checked="" type="checkbox"/>	CAL	California Plant

Analysis, Options, and Filters

Search	Name	Result Status	Filter	Filter (To)
<input checked="" type="checkbox"/>	TPC			
<input checked="" type="checkbox"/>	Prod/Datetime	All		
<input checked="" type="checkbox"/>	Plate Count	All		
<input checked="" type="checkbox"/>	Retest	All		
<input checked="" type="checkbox"/>	Petroleum Lot #	All		

Date Range: **Performed** Beginning: ☒ 1/1/2013 Ending: ☒ 1/31/2013

Location:

Proximity: Zone:

Sampled by User:

Sample Detail

Attn	Analysis	Plant	Location	Plate Count	Lab	Location Description	Sampling Comments	Status	Performed On	Performed By
	TPC	UNK - LINCOLN PLANT	Zone 2-019	1		Overhead Conveyor	Imported from spreadsheet	Closed	1/29/2013 8:...	support
	TPC	UNK - LINCOLN PLANT	Zone 2-019	0		Overhead Conveyor	Imported from spreadsheet	Closed	1/11/2013 8:...	support
	TPC	UNK - LINCOLN PLANT	Zone 2-020	3		Overhead Drop Chutes	Imported from spreadsheet	Closed	1/21/2013 8:...	support
	TPC	UNK - LINCOLN PLANT	Zone 2-020	0		Overhead Drop Chutes	Imported from spreadsheet	Closed	1/2/2013 8:6...	support
	TPC	UNK - LINCOLN PLANT	Zone 2-021	0		Overhead Gates	Imported from spreadsheet	Closed	1/29/2013 8:...	support
	TPC	UNK - LINCOLN PLANT	Zone 2-021	0		Overhead Gates	Imported from spreadsheet	Closed	1/11/2013 8:...	support
	TPC	UNK - LINCOLN PLANT	Zone 2-022	16		Product Rake	Imported from spreadsheet	Closed	1/21/2013 8:...	support
	TPC	UNK - LINCOLN PLANT	Zone 2-022	26		Product Rake	Imported from spreadsheet	Closed	1/2/2013 8:6...	support
	TPC	UNK - LINCOLN PLANT	Zone 2-022	0		Product Rake	Imported from spreadsheet	Closed	1/8/2013 8:6...	support
	TPC	UNK - LINCOLN PLANT	Zone 2-022	4		Product Rake	Imported from spreadsheet	Closed	1/24/2013 8:...	support
	TPC	UNK - LINCOLN PLANT	Zone 2-023	0		Bus Conveyor	Imported from spreadsheet	Closed	1/29/2013 8:...	support
	TPC	UNK - LINCOLN PLANT	Zone 2-023	0		Bus Conveyor	Imported from spreadsheet	Closed	1/11/2013 8:...	support
	TPC	UNK - LINCOLN PLANT	Zone 2-024	0		QA Inspection Area	Imported from spreadsheet	Closed	1/21/2013 8:...	support
	TPC	UNK - LINCOLN PLANT	Zone 2-024	0		QA Inspection Area	Imported from spreadsheet	Closed	1/2/2013 8:6...	support
	TPC	UNK - LINCOLN PLANT	Zone 2-025	7		Cutter	Imported from spreadsheet	Closed	1/26/2013 8:...	support


☒ Show analysis options and values ☒ Show collection detail ☒ Show testing detail

Click on the column title on the table to order the data by a certain criteria. For example, click on the Performed On to see results in chronological order of date performed.

Select and change criteria for the graph using the same panel used for sampling history on the upper portion of the screen.

Hover over points to view details such as date, time, analysis, analysis option, location, etc. This pop-up can be customized.

Hover over a label in the map key to highlight the group of points which belong to that label.



Why isn't an analysis option showing up on my Sampling Graph?

If an option has been checked as **Informational**, it will not appear in the graph.

How to view sample acceptability status over the plant map.

- 1) Click on **Sampling** under the *Main Menu*.
- 2) Enter search criteria and pull up historical data.
- 3) Right-click on one of the rows. A dialog with multiple options pops-up. Click **send all # results to the locations screen for visualization**.




On the left display, locations are highlighted in green and red. Green locations are ones with acceptable counts, and red locations are ones with fail counts.

Highlight a row on the right and the corresponding location will be highlighted automatically on the left display and vice versa.

Order the *Date Sampled* column by clicking on the column title **When**. Click the first row and press the down-arrow on the keyboard to see a play-through of where the samples were taken in a chronological order.

Note: Every row can be ordered by clicking the column title.

Check the box **Color locations based on selected samples only**, and select a group of rows to see the colors only for those selected.

 Clear List Search Samples☒ Color locations based on selected samples only

Part of data not being sent to plant map.

The number of samples transferred to the plant map will depend on the number of items listed in the first page of *Sampling History*. By default, this screen lists 1000 items. If there are more than 1000 samples, there will appear two arrows (one pointing to the left and one to the right) with the text **Next 1000 Samples** on the bottom of the screen.


To increase the limit on number of items to show in this list, click on **Configuration** under the *Main Menu*, then **Settings**.

Click **Sampling** under the *Category* list, and type in 10,000 (more or less) under the *Value* (second column).

Click **Apply**.


How to pull up a sample on the scanner without scanning a barcode.

- 1) Press **Record Analysis**.
- 2) Press **Menu**.
- 3) Press **Manual Input: Tag**.
- 4) Type in the Tag# in the field provided.
- 5) Press **OK**.



How to create a work order on the scanner.

- 1) Press **Work Orders** in the main menu on the scanner.
- 2) Press **Menu**.
- 3) Press **New Work Order**.
- 4) Tap on the first empty box and select an analysis name.
- 5) Tap on the second empty box and select a location name.
- 6) Press **Done** on the scroll pop-up.
- 7) Wait for the newly added sample to appear on the top of the selection area.
- 8) Repeat steps 4, 5, and 6 for each sample.
- 9) Press **Done**.
- 10) Repeat all the steps for each new work order.





What do the result status colors mean on the scanner?

Black = Not Calculated or Error

Red = Unacceptable (out of spec)

Orange = Warning

Green = Acceptable (within spec)

The screenshot shows a mobile application interface titled "Analysis of TPC" with a subtitle "Loading sample 336". At the top are two large buttons: "Done" and "Menu". Below these is a "Sampling Note:" label followed by a text input field. Further down is a "ProdDatetime:" label followed by a text input field containing the word "NOW". Below that is a "Plate Count:" label followed by a text input field containing the number "5", which is highlighted with a green border. Below the plate count is a "Retest:" label followed by a dropdown menu. Further down is a "Petrifilm Lot #:" label followed by a text input field. At the bottom is a "Lab:" label followed by a text input field with a blue information icon (i) to its right.



How to change the date range(reduce clutter) in the work order queue.


There is a menu to the left of the display panel in the *Work Orders* tab.

User is the login user. Make sure to check this box instead of *All* to see work orders specific to you.

Date Range drop-down list shows time periods during which work orders were assigned. For example, selecting “Last Month” will list all work orders created in the previous month. There is an option to set the date range to *Custom*. For this option, users can set a date range by specifying the start and end date. To the right of the date box is a calendar icon for easier date selection.

Select the plant from the drop-down to narrow down the work order list to one plant.

By checking the *Show Closed* box, users can view all completed and pending work orders within the date range set.



Recently created work order is not showing in the *Home Screen* queue.

Click **Refresh** in the *Home Screen* to the left of the work order queue to see the new work order.

Part of data not exported to the .csv file.

The number of samples exported will depend on the number of items listed in the first page of *Sampling History*. By default, this screen lists 1000 items. If there are more than 1000 samples, there will appear two arrows (one pointing to the left and one to the right) with the text **Next 1000 Samples** on the bottom of the screen.

To increase the limit on number of items to show in this list, click on **Configuration** under the *Main Menu*, then **Settings**.

Click **Sampling** under the *Category* list, and type in 10,000 (more or less) under the *Value* (second column).

Click **Apply**.